

EGFR Rabbit anti-Human Polyclonal (aa1010-1029) (FITC) Antibody - LS-C213691 - LSBio	
CatalogID:	LS-C213691
Target:	epidermal growth factor receptor (EGFR)
Synonyms:	EGFR Antibody, EGf receptor Antibody, ERBB Antibody, ERBB1 Antibody, MENA Antibody, Panitumumab Antibody, PIG61 Antibody, Proto-oncogene c-ErbB-1 Antibody, V-erb-b homolog Antibody, EGF-R Antibody, HER1 Antibody
Family / Subfamily:	Protein Kinase / EGF Receptor
Host	EGFR antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	EGFR antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	EGFR antibody was raised against kLH-conjugated synthetic peptide (aa1010-1029) from the C-terminus of human IPO8.
Specificity:	Human EGFR
Epitope:	aa1010-1029
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	