

EGFR Rabbit anti-Human Polyclonal (pSer768) (PE) Antibody - LS-C213714 - LSBio

| | |
|-----------------------------|---|
| CatalogID: | LS-C213714 |
| Target: | epidermal growth factor receptor (EGFR) |
| Synonyms: | EGFR Antibody, EGf receptor Antibody, ERBB Antibody, ERBB1 Antibody, MENA Antibody, Panitumumab Antibody, PIG61 Antibody, Proto-oncogene c-ErbB-1 Antibody, V-erb-b homolog Antibody, EGF-R Antibody, HER1 Antibody |
| Family / Subfamily: | Protein Kinase / EGF Receptor |
| Host | EGFR antibody was produced in Rabbit |
| Clonality: | Polyclonal |
| Isotype: | IgG |
| Conjugations: | Phycoerythrin (PE) |
| Immunogen Species: | EGFR antibody was raised against Human |
| Antigen Type: | Synthetic phospho-peptide |
| Immunogen: | EGFR antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser768 of human EGFR. |
| Specificity: | Human EGFR |
| Epitope: | pSer768 |
| Reactivity: | Human |
| Purification: | Protein A purified |
| Presentation: | PBS, pH 7.2, 0.09% sodium azide. |
| Recommended Storage: | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months. |
| Usage Summary: | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested. |
| Uses: | IHC, ELISA, Dot Blot (Optimal dilution to be determined by the researcher) |
| Size: | 200 µl |
| Requested From: | Japan |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences