

ER Alpha / Estrogen Receptor Rabbit anti-Human Polyclonal (aa130-157) (FITC) Antibody - LS-C213993 - LSBio	
<b>CatalogID:</b>	LS-C213993
<b>Target:</b>	estrogen receptor 1 (ESR1)
<b>Synonyms:</b>	ESR1 Antibody, Era Antibody, ERalpha Antibody, Estrogen receptor Antibody, Estradiol receptor Antibody, ER-alpha Antibody, Estrogen receptor 1 Antibody, NR3A1 Antibody, ER Antibody, ESR Antibody, ESRA Antibody, Estrogen receptor alpha Antibody
<b>Family / Subfamily:</b>	NHR / NR3 Steroid receptor
<b>Host</b>	ESR1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	ER Alpha / Estrogen Receptor antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	ER Alpha / Estrogen Receptor antibody was raised against kLH-conjugated synthetic peptide (aa130-157) from the N-terminal region of human ESR1 isoform1.
<b>Specificity:</b>	Human ER Alpha / Estrogen Receptor
<b>Epitope:</b>	aa130-157
<b>Reactivity:</b>	Human, Mouse
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	