

ER Alpha / Estrogen Receptor Rabbit anti-Human Polyclonal (aa245-273) (FITC) Antibody - LS-C214012 - LSBio	
CatalogID:	LS-C214012
Target:	estrogen receptor 1 (ESR1)
Synonyms:	ESR1 Antibody, Era Antibody, ERalpha Antibody, Estrogen receptor Antibody, Estradiol receptor Antibody, ER-alpha Antibody, Estrogen receptor 1 Antibody, NR3A1 Antibody, ER Antibody, ESR Antibody, ESRA Antibody, Estrogen receptor alpha Antibody
Family / Subfamily:	NHR / NR3 Steroid receptor
Host	ESR1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ER Alpha / Estrogen Receptor antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ER Alpha / Estrogen Receptor antibody was raised against kLH-conjugated synthetic peptide (aa245-273) of human ESR1/ER.
Specificity:	Human ER Alpha / Estrogen Receptor
Epitope:	aa245-273
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/27/2014 © 2014 LifeSpan BioSciences	