

ETS1 / ETS-1 Rabbit anti-Human Polyclonal (aa46-74) (FITC) Antibody - LS-C214021 - LSBio	
<b>CatalogID:</b>	LS-C214021
<b>Target:</b>	v-ets avian erythroblastosis virus E26 oncogene homolog 1 (ETS1)
<b>Synonyms:</b>	ETS1 Antibody, Ets protein Antibody, ETS-1 Antibody, EWSR2 Antibody, p54 Antibody, Protein C-ets-1 Antibody
<b>Family / Subfamily:</b>	Transcription factor / ETS
<b>Host</b>	ETS1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	ETS1 / ETS-1 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	ETS1 / ETS-1 antibody was raised against kLH-conjugated synthetic peptide (aa46-74) from the N-terminal region of human ETS1.
<b>Specificity:</b>	Human ETS1 / ETS-1
<b>Epitope:</b>	aa46-74
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	ICC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	