

EZH2 Rabbit anti-Human Polyclonal (N-Terminus) (FITC) Antibody - LS-C214088 - LSBio

CatalogID:	LS-C214088
Target:	enhancer of zeste homolog 2 (Drosophila) (EZH2)
Synonyms:	EZH2 Antibody, ENX1 Antibody, Enhancer of zeste homolog 2 Antibody, KMT6 Antibody, KMT6A Antibody, Lysine N-methyltransferase 6 Antibody, WVS2 Antibody, Enhancer of zeste Antibody, ENX-1 Antibody
Host	EZH2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	EZH2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	EZH2 antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human EZH2.
Specificity:	Human EZH2.
Epitope:	N-Terminus
Reactivity:	Human, Rabbit
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences