

FAH Rabbit anti-Human Polyclonal (aa14-41) (PE) Antibody - LS-C214125 - LSBio

CatalogID:	LS-C214125
Target:	fumarylacetoacetate hydrolase (fumarylacetoacetase) (FAH)
Synonyms:	FAH Antibody, Beta-diketonase Antibody, Fumarylacetoacetase Antibody, Fumarylacetoacetate hydrolase Antibody, FAA Antibody
Host	FAH antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	FAH antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	FAH antibody was raised against kLH-conjugated synthetic peptide (aa14-41) from the N-terminal region of human FAH.
Specificity:	Human FAH.
Epitope:	aa14-41
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences