

FOSL2 / FRA-2 Rabbit anti-Human Polyclonal (aa95-122) (PE) Antibody - LS-C214310 - LSBio

CatalogID:	LS-C214310
Target:	FOS-like antigen 2 (FOSL2)
Synonyms:	FOSL2 Antibody, FOS-like antigen 2 Antibody, FRA-2 Antibody, Fos-related antigen 2 Antibody, FRA2 Antibody
Host	FOSL2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	FOSL2 / FRA-2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	FOSL2 / FRA-2 antibody was raised against kLH-conjugated synthetic peptide (aa95-122) of human FOSL2.
Specificity:	Human FOSL2 / FRA-2.
Epitope:	aa95-122
Reactivity:	Human, Hamster
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences