

**GCLC Rabbit anti-Human Polyclonal (aa6-34) (PE) Antibody - LS-C214661 - LSBio**

|                             |   |
|-----------------------------|---|
| <b>CatalogID:</b>           | LS-C214661  |
| <b>Target:</b>              | glutamate-cysteine ligase, catalytic subunit (GCLC)   |
| <b>Synonyms:</b>            | GCLC Antibody, GCL Antibody, GCS heavy chain Antibody, Gamma-ECS Antibody, GCS Antibody, GLCL Antibody, GLCLC Antibody                      |
| <b>Host</b>                 | GCLC antibody was produced in Rabbit  |
| <b>Clonality:</b>           | Polyclonal  |
| <b>Isotype:</b>             | IgG   |
| <b>Conjugations:</b>        | Phycoerythrin (PE)  |
| <b>Immunogen Species:</b>   | GCLC antibody was raised against Human  |
| <b>Antigen Type:</b>        | Synthetic peptide   |
| <b>Immunogen:</b>           | GCLC antibody was raised against kLH-conjugated synthetic peptide (aa6-34) from the N-terminal region of human GCLC.                        |
| <b>Specificity:</b>         | Human GCLC.   |
| <b>Epitope:</b>             | aa6-34  |
| <b>Reactivity:</b>          | Human   |
| <b>Purification:</b>        | Protein A purified  |
| <b>Presentation:</b>        | PBS, pH 7.2, 0.09% sodium azide.  |
| <b>Recommended Storage:</b> | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months. |
| <b>Usage Summary:</b>       | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.                  |
| <b>Uses:</b>                | IHC, ICC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)   |
| <b>Size:</b>                | 200 µl  |
| <b>Requested From:</b>      | Japan   |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences