

GCLM Rabbit anti-Human Polyclonal (aa245-274) (PE) Antibody - LS-C214664 - LSBio

CatalogID:	LS-C214664
Target:	glutamate-cysteine ligase, modifier subunit (GCLM)
Synonyms:	GCLM Antibody, Gamma-ECS regulatory subunit Antibody, GLCLR Antibody, GCS light chain Antibody, GSC light chain Antibody
Host	GCLM antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	GCLM antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	GCLM antibody was raised against kLH-conjugated synthetic peptide (aa245-274) from the C-terminus of human GCLM.
Specificity:	Human GCLM.
Epitope:	aa245-274
Reactivity:	Human, Mouse
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences