

GGT1 / GGT Rabbit anti-Mouse Polyclonal (aa12-39) (PE) Antibody - LS-C214695 - LSBio	
CatalogID:	LS-C214695
Target:	gamma-glutamyltransferase 1 (GGT1)
Synonyms:	GGT1 Antibody, CD224 antigen Antibody, CD224 Antibody, Ec 2.3.2.2 Antibody, Glutamyl transpeptidase Antibody, GTG Antibody, Gamma-glutamyltransferase Antibody, Glutathione hydrolase 1 Antibody, GLUTATHIONURIA Antibody, Gamma-glutamyl transpeptidase Antibody, Gamma-glutamyltransferase 1 Antibody, Gamma-glutamyltranspeptidase 1 Antibody, GGT Antibody, Leukotriene-C4 hydrolase Antibody, D22S672 Antibody, D22S732 Antibody, Gamma-glutamyltransferase 2 Antibody, GGT 1 Antibody
Family / Subfamily:	Protease / Threonine T3
Host	GGT1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	GGT1 / GGT antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	GGT1 / GGT antibody was raised against kLH-conjugated synthetic peptide (aa12-39) from the N-terminal region of rat Ggt1.
Specificity:	Mouse GGT1 / GGT.
Epitope:	aa12-39
Reactivity:	Mouse, Rabbit
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	