

GLMN Rabbit anti-Human Polyclonal (aa505-533) (PE) Antibody - LS-C214724 - LSBio

CatalogID:	LS-C214724
Target:	glomulin, FKBP associated protein (GLMN)
Synonyms:	GLMN Antibody, FAP68 Antibody, GLML Antibody, GVM Antibody, FKBP-associated protein Antibody, Glomulin Antibody, VMGLOM Antibody, FAP48 Antibody, FKBPAP Antibody
Host	GLMN antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	GLMN antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	GLMN antibody was raised against kLH-conjugated synthetic peptide (aa505-533) from the C-terminus of human GLMN.
Specificity:	Human GLMN.
Epitope:	aa505-533
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences