

CD3E Hamster anti-Mouse Monoclonal (PE, Cy5) (145-2C11) Antibody - LS-C21527 - LSBio

CatalogID:	LS-C21527
Target:	CD3e molecule, epsilon (CD3-TCR complex)
Synonyms:	CD3E Antibody, CD3e antigen Antibody, CD3-epsilon Antibody, T-cell receptor epsilon chain Antibody, TCRE Antibody, T-cell receptors epsilon chain Antibody, T3E Antibody, CD3 Epsilon Antibody
Host	CD3E antibody was produced in Hamster
Clonality:	Monoclonal
Isotype:	IgG,k
Clone Name:	145-2C11
Conjugations:	Phycoerythrin (PE), Cy5
Immunogen Species:	CD3E antibody was raised against Mouse
Immunogen:	CD3E antibody was raised against h-2Kb-specific mouse cytotoxic T lymphocyte clone BM10-37.
Specificity:	Mouse CD3e (Mr 25kD). Recognizes an epitope on the 25kD epsilon chain of the CD3/TCR complex.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, pH 7.4, sodium azide, a stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Identification and enumeration of CD3+ cells by flow cytometry. Immunoprecipitation. Western blotting. In vitro depletion of CD3+ cells. In vitro activation of T cells. Immunohistochemistry. Flow Cytometry: 1 ug/1E6 cells.
Uses:	IHC, Western blot, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences