

**CD3E Hamster anti-Mouse Monoclonal (PE, Cy5) Antibody - LS-C21535 - LSBio**

<b>CatalogID:</b>	LS-C21535
<b>Target:</b>	CD3e molecule, epsilon (CD3-TCR complex)
<b>Synonyms:</b>	CD3E Antibody, CD3e antigen Antibody, CD3-epsilon Antibody, T-cell receptor epsilon chain Antibody, TCRE Antibody, T-cell receptors epsilon chain Antibody, T3E Antibody, CD3 Epsilon Antibody
<b>Host</b>	CD3E antibody was produced in Hamster
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG,k
<b>Conjugations:</b>	Phycoerythrin (PE), Cy5
<b>Immunogen Species:</b>	CD3E antibody was raised against Mouse
<b>Immunogen:</b>	CD3E antibody was raised against t cell receptor complexes derived from C6VL-BS thymoma cells (C57BL/Ka).
<b>Specificity:</b>	Mouse CD3e, Mr 25kD. Recognizes an epitope on the 25kD epsilon chain of the CD3/TCR complex of mouse strains tested.
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, pH 7.4, sodium azide, a stabilizing agent.
<b>Recommended Storage:</b>	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
<b>Usage Summary:</b>	Identification and enumeration of CD3+ cells by flow cytometry. Immunoprecipitation. In vitro activation of T cells. Flow Cytometry: 1 ug/1E6 cells.
<b>Uses:</b>	Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences