

CD3E Rat anti-Mouse Monoclonal (PE) Antibody - LS-C21541 - LSBio

CatalogID:	LS-C21541
Target:	CD3e molecule, epsilon (CD3-TCR complex)
Synonyms:	CD3E Antibody, CD3e antigen Antibody, CD3-epsilon Antibody, T-cell receptor epsilon chain Antibody, TCRE Antibody, T-cell receptors epsilon chain Antibody, T3E Antibody, CD3 Epsilon Antibody
Host	CD3E antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2b,k
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD3E antibody was raised against Mouse
Specificity:	Mouse CD3e (Mr 25kD). Recognizes an epitope on the 25kD epsilon chain of the CD3/TCR complex.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, pH 7.4, sodium azide, a stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Suitable for use in in vitro depletion of CD3+ cells, in vitro activation of T cells (accessory cells must be provided), Immunoprecipitation, Flow Cytometry and Immunohistochemistry. Flow Cytometry: 1 ug/1E6 cells.
Uses:	IHC, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Concentration:	0.1 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences