

CD3E Mouse anti-Pig Monoclonal (FITC) Antibody - LS-C21551 - LSBio	
CatalogID:	LS-C21551
Target:	CD3e molecule, epsilon (CD3-TCR complex)
Synonyms:	CD3E Antibody, CD3e antigen Antibody, CD3-epsilon Antibody, T-cell receptor epsilon chain Antibody, TCRE Antibody, T-cell receptors epsilon chain Antibody, T3E Antibody, CD3 Epsilon Antibody
Host	CD3E antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD3E antibody was raised against Pig
Specificity:	Recognizes porcine CD3e, Mr 20kD.
Reactivity:	Pig
Purification:	Purified
Presentation:	PBS, sodium azide.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. Immunohistochemistry (acetone-fixed, frozen tissue sections). Immunoprecipitation. In vitro activation of T cells. Flow Cytometry: 1 ug/1E6 cells.
Uses:	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	