

**CD3E Mouse anti-Pig Monoclonal (PE, Cy5) (PPT3) Antibody - LS-C21554 - LSBio**

<b>CatalogID:</b>	LS-C21554
<b>Target:</b>	CD3e molecule, epsilon (CD3-TCR complex)
<b>Synonyms:</b>	CD3E Antibody, CD3e antigen Antibody, CD3-epsilon Antibody, T-cell receptor epsilon chain Antibody, TCRE Antibody, T-cell receptors epsilon chain Antibody, T3E Antibody, CD3 Epsilon Antibody
<b>Host</b>	CD3E antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1,k
<b>Clone Name:</b>	PPT3
<b>Conjugations:</b>	Phycoerythrin (PE), Cy5
<b>Immunogen Species:</b>	CD3E antibody was raised against Pig
<b>Specificity:</b>	Recognizes porcine CD3e, Mr 20kD.
<b>Reactivity:</b>	Pig
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, 0.1% sodium azide, stabilizing agent.
<b>Recommended Storage:</b>	Store at +4°C. DO NOT FREEZE.
<b>Usage Summary:</b>	Flow cytometry. Immunohistochemistry (acetone-fixed, frozen tissue sections). Immunoprecipitation. In vitro activation of T cells. Flow Cytometry: 0.2 ug/1E6 cells.
<b>Uses:</b>	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Concentration:</b>	0.1 mg/ml
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences