

CD3E Mouse anti-Pig Monoclonal (FITC) Antibody - LS-C21556 - LSBio	
<b>CatalogID:</b>	LS-C21556
<b>Target:</b>	CD3e molecule, epsilon (CD3-TCR complex)
<b>Synonyms:</b>	CD3E Antibody, CD3e antigen Antibody, CD3-epsilon Antibody, T-cell receptor epsilon chain Antibody, TCRE Antibody, T-cell receptors epsilon chain Antibody, T3E Antibody, CD3 Epsilon Antibody
<b>Host</b>	CD3E antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2b,k
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	CD3E antibody was raised against Pig
<b>Specificity:</b>	Recognizes porcine CD3e, Mr 20kD.
<b>Reactivity:</b>	Pig
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, sodium azide.
<b>Recommended Storage:</b>	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
<b>Usage Summary:</b>	Flow cytometry. Immunoprecipitation. Flow Cytometry: 1 ug/1E6 cells.
<b>Uses:</b>	Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	500 µg
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b> <b>Not for resale without prior written consent from LifeSpan BioSciences, Inc.</b> Created on 9/27/2014 © 2014 LifeSpan BioSciences	