

MAP1LC3B / LC3B Rabbit anti-Human Polyclonal (Thr93/Tyr99) (FITC) Antibody - LS-C215689 - LSBio

<b>CatalogID:</b>	LS-C215689
<b>Target:</b>	microtubule-associated protein 1 light chain 3 beta (MAP1LC3B)
<b>Synonyms:</b>	MAP1LC3B Antibody, ATG8F Antibody, LC3B Antibody, MAP1A/1BLC3 Antibody, MAP1A/MAP1B LC3 B Antibody, MAP1A/MAP1B light chain 3 B Antibody, MAP1ALC3 Antibody, MAP1LC3B-a Antibody
<b>Host</b>	MAP1LC3B antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	MAP1LC3B / LC3B antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	MAP1LC3B / LC3B antibody was raised against kLH-conjugated synthetic peptide from around Thr93/Tyr99 of human APG8b.
<b>Specificity:</b>	Human MAP1LC3B / LC3B
<b>Epitope:</b>	Thr93/Tyr99
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences