

MPO / Myeloperoxidase Rabbit anti-Human Polyclonal (aa419-448) (FITC) Antibody - LS-C215883 - LSBio	
CatalogID:	LS-C215883
Target:	myeloperoxidase (MPO)
Synonyms:	MPO Antibody, Myeloperoxidase precursor Antibody, Myeloperoxidase Antibody
Host	MPO antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	MPO / Myeloperoxidase antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MPO / Myeloperoxidase antibody was raised against kLH-conjugated synthetic peptide (aa419-448) of human Myeloperoxidase.
Specificity:	Human MPO / Myeloperoxidase
Epitope:	aa419-448
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	