

PIN1 Rabbit anti-Human Polyclonal (aa36-66) (FITC) Antibody - LS-C216782 - LSBio

CatalogID:	LS-C216782
Target:	peptidylprolyl cis/trans isomerase, NIMA-interacting 1 (PIN1)
Synonyms:	PIN1 Antibody, DOD Antibody, PPlase Pin1 Antibody, Rotamase Pin1 Antibody
Host	PIN1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	PIN1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PIN1 antibody was raised against kLH-conjugated synthetic peptide (aa36-66) from internal region of human PIN1.
Specificity:	Human PIN1.
Epitope:	aa36-66
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences