

**PNCK Rabbit anti-Human Polyclonal (aa314-343) (FITC) Antibody - LS-C216872 - LSBio**

<b>CatalogID:</b>	LS-C216872
<b>Target:</b>	pregnancy up-regulated nonubiquitous CaM kinase (PNCK)
<b>Synonyms:</b>	PNCK Antibody, BSTK3 Antibody, CaM kinase I beta Antibody, CaM kinase IB Antibody, CaM kinase Ibeta1 Antibody, CaM-KI beta Antibody, CaMK1b Antibody, CaMKI-beta Antibody, CaM kinase Ibeta2 Antibody, CaMKIbeta2 Antibody
<b>Family / Subfamily:</b>	Protein Kinase / CAMK1
<b>Host</b>	PNCK antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	PNCK antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	PNCK antibody was raised against kLH-conjugated synthetic peptide (aa314-343) from the C-terminus of human PNCK.
<b>Specificity:</b>	Human PNCK
<b>Epitope:</b>	aa314-343
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences