

**CD28 Hamster anti-Mouse Monoclonal (FITC) Antibody - LS-C21701 - LSBio**

<b>CatalogID:</b>	LS-C21701
<b>Target:</b>	CD28 molecule
<b>Synonyms:</b>	CD28 Antibody, CD28 antigen Antibody, CD28 antigen (Tp44) Antibody, Tp44 Antibody, CD28 molecule Antibody
<b>Host</b>	CD28 antibody was produced in Hamster
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	CD28 antibody was raised against Mouse
<b>Specificity:</b>	Recognizes murine CD28.
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, sodium azide.
<b>Recommended Storage:</b>	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
<b>Usage Summary:</b>	Flow cytometry. Immunoprecipitation. In vitro costimulation of T cells and NK cells. Flow Cytometry: 3 ug/1E6 cells.
<b>Uses:</b>	Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	500 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences