

RNPEP Rabbit anti-Human Polyclonal (aa380-411) (FITC) Antibody - LS-C217238 - LSBio

CatalogID:	LS-C217238
Target:	Arginyl Aminopeptidase (aminopeptidase B) (RNPEP)
Synonyms:	RNPEP Antibody, APB Antibody, Arginine aminopeptidase Antibody, Arginyl aminopeptidase Antibody, AP-B Antibody, Aminopeptidase B Antibody
Family / Subfamily:	Exopeptidase / Metallopeptidase M1
Host	RNPEP antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	RNPEP antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	RNPEP antibody was raised against kLH-conjugated synthetic peptide (aa380-411) of human RNPEP.
Specificity:	Human RNPEP
Epitope:	aa380-411
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences