

CD28 Mouse anti-Chicken Monoclonal (PE) Antibody - LS-C21726 - LSBio

CatalogID:	LS-C21726
Target:	CD28 molecule
Synonyms:	CD28 Antibody, CD28 antigen Antibody, CD28 antigen (Tp44) Antibody, Tp44 Antibody, CD28 molecule Antibody
Host	CD28 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD28 antibody was raised against Chicken
Specificity:	Recognizes chicken CD28, a monomeric molecule of Mr 40-44kD with 50% amino acid sequence identity with mammalian CD28. Also reacts with turkey CD28 as demonstrated by flow cytometry and immunoprecipitation/Western blot.
Reactivity:	Chicken
Purification:	Purified
Presentation:	PBS, pH 7.4, 0.1% sodium azide, a stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. Immunoprecipitation. In vitro T-cell costimulation. Flow Cytometry: 0.2 ug/1E6 cells.
Uses:	Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences