

SIRT3 / Sirtuin 3 Rabbit anti-Mouse Polyclonal (aa229-259) (PE) Antibody - LS-C217513 - LSBio	
CatalogID:	LS-C217513
Target:	sirtuin 3 (SIRT3)
Synonyms:	SIRT3 Antibody, HSIRT3 Antibody, Sir2-like 3 Antibody, SIR2-like protein 3 Antibody, SIR2L3 Antibody, Sirtuin 3 Antibody, Sirtuin type 3 Antibody
Family / Subfamily:	ADP-Ribosyltransferase / not assigned-ADP-Ribosyltransferase
Host	SIRT3 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	SIRT3 / Sirtuin 3 antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	SIRT3 / Sirtuin 3 antibody was raised against kLH-conjugated synthetic peptide (aa229-259) from the C-terminus of mouse Sirt3.
Specificity:	Mouse SIRT3
Epitope:	aa229-259
Reactivity:	Mouse, Human, Hamster
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	