

CD33 Mouse anti-Human Monoclonal (PE) Antibody - LS-C21761 - LSBio

CatalogID:	LS-C21761
Target:	CD33 molecule
Synonyms:	CD33 Antibody, CD33 antigen Antibody, CD33 molecule Antibody, gp67 Antibody, SIGLEC-3 Antibody, p67 Antibody, CD33 antigen (gp67) Antibody, SIGLEC3 Antibody
Family / Subfamily:	Immunoglobulin / not assigned-Immunoglobulin
Host	CD33 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD33 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD33 antibody was raised against human AML cells.
Specificity:	Recognizes the human CD33 cell surface glycoprotein, a 67kD molecule expressed by monocytes, granulocyte precursors, myeloid progenitor cells and myeloid leukemias.
Reactivity:	Human
Purification:	Protein G purified
Presentation:	Lyophilized from PBS, pH 7.2, 0.1% sodium azide, 1% BSA.
Usage Summary:	Flow Cytometry: Neat.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences