

CD33 Mouse anti-Human Monoclonal (PE, Cy5) Antibody - LS-C21773 - LSBio

CatalogID:	LS-C21773
Target:	CD33 molecule
Synonyms:	CD33 Antibody, CD33 antigen Antibody, CD33 molecule Antibody, gp67 Antibody, SIGLEC-3 Antibody, p67 Antibody, CD33 antigen (gp67) Antibody, SIGLEC3 Antibody
Family / Subfamily:	Immunoglobulin / not assigned-Immunoglobulin
Host	CD33 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Phycoerythrin (PE), Cy5
Immunogen Species:	CD33 antibody was raised against Human
Specificity:	Recognizes human CD33.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. Immunoprecipitation. Western blotting. Immunohistochemistry (frozen sections). Flow Cytometry: 10 ul/1E6 cells.
Uses:	IHC - Frozen, Western blot, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences