

TANK Rabbit anti-Human Polyclonal (aa95-122) (PE) Antibody - LS-C217838 - LSBio

CatalogID:	LS-C217838
Target:	TRAF family member-associated NFKB activator (TANK)
Synonyms:	TANK Antibody, I-TRAF Antibody, TRAF-interacting protein Antibody, ITRAF Antibody
Host	TANK antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	TANK antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	TANK antibody was raised against kLH-conjugated synthetic peptide (aa95-122) of human TANK.
Specificity:	Human TANK
Epitope:	aa95-122
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences