

**CD45R Rat anti-Mouse Monoclonal (FITC) Antibody - LS-C21790 - LSBio**

<b>CatalogID:</b>	LS-C21790
<b>Target:</b>	CD45R
<b>Host</b>	CD45R antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	CD45R antibody was raised against Mouse
<b>Specificity:</b>	Reacts with mouse CD45R (Ly-5, B220) expressed on pre-B cells and B cells.
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	Purified IgG.
<b>Usage Summary:</b>	Flow Cytometry: Neat-1:10.
<b>Uses:</b>	Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences