

CD45R Rat anti-Mouse Monoclonal (PE) Antibody - LS-C21796 - LSBio	
CatalogID:	LS-C21796
Target:	CD45R
Host	CD45R antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a,k
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD45R antibody was raised against Mouse
Immunogen:	CD45R antibody was raised against abelson murine leukemia virus-induced pre-B tumor cells.
Specificity:	Recognizes B220 isoform of CD45, Mr 240kD.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. Immunohistochemistry (frozen sections). Immunoprecipitation. Flow Cytometry: 0.2 ug/1E6 cells.
Uses:	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/28/2014 © 2014 LifeSpan BioSciences	