

CD45R Rat anti-Mouse Monoclonal (PE, Cy5) Antibody - LS-C21798 - LSBio	
<b>CatalogID:</b>	LS-C21798
<b>Target:</b>	CD45R
<b>Host</b>	CD45R antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a,k
<b>Conjugations:</b>	Phycoerythrin (PE), Cy5
<b>Immunogen Species:</b>	CD45R antibody was raised against Mouse
<b>Immunogen:</b>	CD45R antibody was raised against abelson murine leukemia virus-induced pre-B tumor cells.
<b>Specificity:</b>	Recognizes B220 isoform of CD45, Mr 240kD.
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, 0.1% sodium azide, stabilizing agent.
<b>Recommended Storage:</b>	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
<b>Usage Summary:</b>	Flow cytometry. Immunohistochemistry (frozen sections). Immunoprecipitation. Flow Cytometry: 0.2 ug/1E6 cells.
<b>Uses:</b>	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b> Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/28/2014 © 2014 LifeSpan BioSciences	