

CD45R Rat anti-Mouse Monoclonal Antibody - LS-C21802 - LSBio

CatalogID:	LS-C21802
Target:	CD45R
Host	CD45R antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a,k
Immunogen Species:	CD45R antibody was raised against Mouse
Immunogen:	CD45R antibody was raised against abelson murine leukemia virus-induced pre-B tumor cells.
Specificity:	Recognizes B220 isoform of CD45, Mr 240kD.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS. No preservative added.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Suitable for use in Flow Cytometry, Immunoprecipitation and Immunohistochemistry. Immunohistochemistry: frozen.
Uses:	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	500 µg
Concentration:	0.5 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences