

CD45RA Mouse anti-Human Monoclonal (PE, Cy5) Antibody - LS-C21818 - LSBio

<b>CatalogID:</b>	LS-C21818
<b>Target:</b>	CD45RA
<b>Host</b>	CD45RA antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Conjugations:</b>	Phycoerythrin (PE), Cy5
<b>Immunogen Species:</b>	CD45RA antibody was raised against Human
<b>Specificity:</b>	Recognizes human CD45RA.
<b>Reactivity:</b>	Human
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, 0.1% sodium azide, stabilizing agent.
<b>Recommended Storage:</b>	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
<b>Usage Summary:</b>	Flow cytometry. Immunohistochemistry (frozen sections and paraffin sections). Immunoprecipitation. Flow Cytometry: 10 ul/1E6 cells.
<b>Uses:</b>	IHC - Paraffin, IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 tst
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences