

ULK1 Rabbit anti-Human Polyclonal (pSer556) (FITC) Antibody - LS-C218316 - LSBio

CatalogID:	LS-C218316
Target:	unc-51 like autophagy activating kinase 1 (ULK1)
Synonyms:	ULK1 Antibody, ATG1A Antibody, ATG1 Antibody, HATG1 Antibody, KIAA0722 Antibody, Unc-51-like kinase 1 Antibody, Unc51.1 Antibody, UNC51 Antibody
Family / Subfamily:	Protein Kinase / ULK
Host	ULK1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ULK1 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	ULK1 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser556 of human ULK1.
Specificity:	Human ULK1.
Epitope:	pSer556
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences