

CD45RA Mouse anti-Human Monoclonal Antibody - LS-C21834 - LSBio	
CatalogID:	LS-C21834
Target:	CD45RA
Host	CD45RA antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Immunogen Species:	CD45RA antibody was raised against Human
Immunogen:	CD45RA antibody was raised against human CD45RA.
Specificity:	Recognizes human CD45RA antigen.
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.01% sodium azide.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Suitable for use in Flow Cytometry, Immunohistochemistry, Fluorescent Microscopy and Western Blot. Immunohistochemistry: Frozen sections.
Uses:	IHC - Frozen, Western blot, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	