

CD45RA Mouse anti-Human Monoclonal Antibody - LS-C21835 - LSBio	
CatalogID:	LS-C21835
Target:	CD45RA
Host	CD45RA antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Immunogen Species:	CD45RA antibody was raised against Human
Immunogen:	CD45RA antibody was raised against cD45RA.
Specificity:	Recognizes human CD45RA (B-cells, T-cells, NK-cells).
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.01% sodium azide.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Suitable for use in Flow Cytometry.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	