

CD45RB Mouse anti-Human Monoclonal Antibody - LS-C21859 - LSBio	
<b>CatalogID:</b>	LS-C21859
<b>Target:</b>	CD45RB
<b>Host</b>	CD45RB antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Immunogen Species:</b>	CD45RB antibody was raised against Human
<b>Specificity:</b>	Recognizes the isoform of the leucocyte common antigen encoded by exon B and stains B cells, most T cells and monocytes. Clone 3H1434 immunoprecipitates bands of 190, 200 and 205kD. Species cross-reactivity: Rhesus monkey.
<b>Reactivity:</b>	Human, Monkey
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2, 0.1% sodium azide.
<b>Usage Summary:</b>	Immunohistochemistry: Frozen, Paraffin. Immunoprecipitation. Flow Cytometry: 1:50-1:100.
<b>Uses:</b>	IHC - Paraffin, IHC - Frozen, Immunoprecipitation, Flow Cytometry (1:50 - 1:100) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µg
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	