

AKT2 Rabbit anti-Human Polyclonal (pSer474) (PE) Antibody - LS-C218598 - LSBio

<b>CatalogID:</b>	LS-C218598
<b>Target:</b>	v-akt murine thymoma viral oncogene homolog 2 (AKT2)
<b>Synonyms:</b>	AKT2 Antibody, PKB beta Antibody, PKBBETA Antibody, PRKBB Antibody, Protein kinase Akt-2 Antibody, Protein kinase B beta Antibody, RAC-PK-beta Antibody, Rac protein kinase beta Antibody, RAC-BETA Antibody, HIHGHH Antibody, PKBB Antibody
<b>Family / Subfamily:</b>	Protein Kinase / RAC/AKT
<b>Host</b>	AKT2 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	AKT2 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic phospho-peptide
<b>Immunogen:</b>	AKT2 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser474 of human AKT2.
<b>Specificity:</b>	Human AKT2.
<b>Epitope:</b>	pSer474
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Western blot, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences