

AKT2 Rabbit anti-Human Polyclonal (N-Terminus) (FITC) Antibody - LS-C218602 - LSBio

CatalogID:	LS-C218602
Target:	v-akt murine thymoma viral oncogene homolog 2 (AKT2)
Synonyms:	AKT2 Antibody, PKB beta Antibody, PKBBETA Antibody, PRKBB Antibody, Protein kinase Akt-2 Antibody, Protein kinase B beta Antibody, RAC-PK-beta Antibody, Rac protein kinase beta Antibody, RAC-BETA Antibody, HIHGHH Antibody, PKBB Antibody
Family / Subfamily:	Protein Kinase / RAC/AKT
Host	AKT2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	AKT2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	AKT2 antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human AKT2.
Specificity:	Human AKT2.
Epitope:	N-Terminus
Reactivity:	Human, Mouse
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences