

APP / Beta Amyloid Rabbit anti-Human Polyclonal (pSer730) (FITC) Antibody - LS-C218639 - LSBio

<b>CatalogID:</b>	LS-C218639
<b>Target:</b>	amyloid beta (A4) precursor protein (APP)
<b>Synonyms:</b>	APP Antibody, ABPP Antibody, A4 Antibody, ABETA Antibody, AD1 Antibody, Alzheimer disease Antibody, Amyloid beta Antibody, Beta-amyloid peptide Antibody, APPI Antibody, Beta amyloid Antibody, CVAP Antibody, AAA Antibody, PN2 Antibody, PreA4 Antibody, Amyloid beta A4 protein Antibody, Peptidase nexin -II Antibody, CTFgamma Antibody, PN-II Antibody, Protease nexin-II Antibody
<b>Family / Subfamily:</b>	Amyloid precursor protein / Not assigned-Other
<b>Host</b>	APP antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	APP / Beta Amyloid antibody was raised against Human
<b>Antigen Type:</b>	Synthetic phospho-peptide
<b>Immunogen:</b>	APP / Beta Amyloid antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser730 of human APP.
<b>Specificity:</b>	Human APP / Beta Amyloid
<b>Epitope:</b>	pSer730
<b>Reactivity:</b>	Human
<b>Purification:</b>	Immunoaffinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences