

APP / Beta Amyloid Rabbit anti-Human Polyclonal (C-Terminus) (FITC) Antibody - LS-C218650 - LSBio

CatalogID:	LS-C218650
Target:	amyloid beta (A4) precursor protein (APP)
Synonyms:	APP Antibody, ABPP Antibody, A4 Antibody, ABETA Antibody, AD1 Antibody, Alzheimer disease Antibody, Amyloid beta Antibody, Beta-amyloid peptide Antibody, APPI Antibody, Beta amyloid Antibody, CVAP Antibody, AAA Antibody, PN2 Antibody, PreA4 Antibody, Amyloid beta A4 protein Antibody, Peptidase nexin -II Antibody, CTFgamma Antibody, PN-II Antibody, Protease nexin-II Antibody
Family / Subfamily:	Amyloid precursor protein / Not assigned-Other
Host	APP antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	APP / Beta Amyloid antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	APP / Beta Amyloid antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human APP.
Specificity:	Human APP / Beta Amyloid
Epitope:	C-Terminus
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences