

CDH1 / E Cadherin Rabbit anti-Human Polyclonal (C-Terminus) (Azide-free) (HRP) Antibody - LS-C218710 - LSBio	
CatalogID:	LS-C218710
Target:	cadherin 1, type 1, E-cadherin (epithelial) (CDH1)
Synonyms:	CDH1 Antibody, CD324 Antibody, CDHE Antibody, Cell-CAM 120/80 Antibody, E Cadherin Antibody, E-Cadherin Antibody, ECAD Antibody, Epithelial cadherin Antibody, Arc-1 Antibody, UVO Antibody, Cadherin-1 Antibody, CAM 120/80 Antibody, CD324 antigen Antibody, LCAM Antibody, Uvomorulin Antibody
Family / Subfamily:	Cadherin / not assigned-Cadherin
Host	CDH1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Horseradish Peroxidase (HRP)
Modifications:	Azide-free
Immunogen Species:	CDH1 / E Cadherin antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CDH1 / E Cadherin antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human CDH1.
Specificity:	Human CDH1 / E Cadherin
Epitope:	C-Terminus
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2. No preservative added (azide free)
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/28/2014 © 2014 LifeSpan BioSciences	