

CD45RC Rat anti-Mouse Monoclonal Antibody - LS-C21886 - LSBio	
CatalogID:	LS-C21886
Target:	CD45RC
Host	CD45RC antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a,k
Immunogen Species:	CD45RC antibody was raised against Mouse
Specificity:	Recognizes an exon C-dependent epitope of CD45 glycoprotein.
Reactivity:	Mouse
Purification:	Purified
Presentation:	100 mM borate buffered saline, pH 8.0. No preservatives or amine-containing buffer salts added.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Identification and enumeration of CD45RC+ cells by flow cytometry. Immunoprecipitation.
Uses:	Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	