

CD45RC Rat anti-Mouse Monoclonal (Biotin) Antibody - LS-C21888 - LSBio

CatalogID:	LS-C21888
Target:	CD45RC
Host	CD45RC antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a,k
Conjugations:	Biotin
Immunogen Species:	CD45RC antibody was raised against Mouse
Specificity:	Recognizes an exon C-dependent epitope of CD45 glycoprotein.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, sodium azide.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Identification and enumeration of CD45RC+ cells by flow cytometry. Immunoprecipitation. Flow Cytometry: 1 ug/1E6 cells.
Uses:	Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences