

CD45RC Rat anti-Mouse Monoclonal (Biotin) Antibody - LS-C21893 - LSBio	
<b>CatalogID:</b>	LS-C21893
<b>Target:</b>	CD45RC
<b>Host</b>	CD45RC antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgM,k
<b>Conjugations:</b>	Biotin
<b>Immunogen Species:</b>	CD45RC antibody was raised against Mouse
<b>Specificity:</b>	Recognizes an exon C-dependent epitope of CD45 glycoprotein.
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, sodium azide.
<b>Recommended Storage:</b>	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
<b>Usage Summary:</b>	Flow Cytometry: 1 ug/1E6 cells.
<b>Uses:</b>	Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	500 µg
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	