

EGFR Rabbit anti-Human Polyclonal (pTyr1016) (PE) Antibody - LS-C219046 - LSBio

CatalogID:	LS-C219046
Target:	epidermal growth factor receptor (EGFR)
Synonyms:	EGFR Antibody, EGf receptor Antibody, ERBB Antibody, ERBB1 Antibody, MENA Antibody, Panitumumab Antibody, PIG61 Antibody, Proto-oncogene c-ErbB-1 Antibody, V-erb-b homolog Antibody, EGF-R Antibody, HER1 Antibody
Family / Subfamily:	Protein Kinase / EGF Receptor
Host	EGFR antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	EGFR antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	EGFR antibody was raised against kLH-conjugated synthetic phosphopeptide around Tyr1016 of human EGFR.
Specificity:	Human EGFR
Epitope:	pTyr1016
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences