

CD95 / FAS Rabbit anti-Human Polyclonal (pTyr291) (FITC) Antibody - LS-C219228 - LSBio

CatalogID:	LS-C219228
Target:	Fas cell surface death receptor
Synonyms:	FAS Antibody, ALPS1A Antibody, APO-1 Antibody, APT1 Antibody, Apo-1 antigen Antibody, APO-1 cell surface antigen Antibody, Apo11 antigen Antibody, Apoptosis antigen 1 Antibody, CD95 antigen Antibody, Delta Fas/APO-1/CD95 Antibody, Fas AMA Antibody, FAS1 Antibody, FasI receptors Antibody, FASLG receptor Antibody, FASTM Antibody, FAS 827dupA Antibody, FasI receptor Antibody, TNFRSF6 Antibody, CD95 Antibody
Family / Subfamily:	TNF Receptor / not assigned-TNF Receptor
Host	FAS antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD95 / FAS antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	CD95 / FAS antibody was raised against kLH-conjugated synthetic phosphopeptide around Tyr291 of human FAS.
Specificity:	Human CD95 / FAS
Epitope:	pTyr291
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences