

FDPS Rabbit anti-Human Polyclonal (Internal) (FITC) Antibody - LS-C219276 - LSBio

CatalogID:	LS-C219276
Target:	farnesyl diphosphate synthase (FDPS)
Synonyms:	FDPS Antibody, Dimethylallyltranstransferase Antibody, Farnesyl diphosphate synthase Antibody, FPP synthetase Antibody, FPS Antibody, KIAA1293 Antibody, FPPS Antibody, Geranyltranstransferase Antibody, FPP synthase Antibody
Host	FDPS antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	FDPS antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	FDPS antibody was raised against kLH-conjugated synthetic peptide from internal region of human FDPS.
Specificity:	Human FDPS.
Epitope:	Internal
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences