

FPGS Rabbit anti-Human Polyclonal (Internal) (PE) Antibody - LS-C219364 - LSBio

CatalogID:	LS-C219364
Target:	folylpolyglutamate synthase (FPGS)
Synonyms:	FPGS Antibody, Folylpolyglutamate synthase Antibody, Tetrahydrofolate synthase Antibody
Host	FPGS antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	FPGS antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	FPGS antibody was raised against kLH-conjugated synthetic peptide from internal region of human FPGS.
Specificity:	Human FPGS
Epitope:	Internal
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences